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1.Introduction

Thank you for purchasing this IR Thermometer. This thermometer is a non-contact infrared measuring instrument designed for temperature evaluation of forehead and household appliances. Features include a 4 digits backlit LCD, scan/hold function and auto power off (20 seconds.) To measure a temperature, point the unit at the object, pull the trigger to scan and release the trigger to hold. Make sure the target area is larger than the unit's spot size.


2.Safety Information

Read the following safety information carefully before attempting to operate or service the meter. Only qualified personnel should perform repairs or servicing not covered in this manual.

2-1 Cautions & Warnings


As with any thermometer, proper technique is crucial to getting accurate temperature readings. Please read this manual thoroughly and carefully before using.

- Do not submerge the unit in water.
- Do not disassemble the thermometer .
- This product is not designed for use in medical evaluations.
- Avoid dropping the thermometer from a height or strongly hitting it with a hard object.
- Always store the thermometer in a cool and dry place: temperatures between -10°C (14°F) and 60°C (140°F); relative humidity less than 95% avoid sunlight.
- The product can only be used to measure forehead temperature simply for reference. If the reading is high (for example, above 38°C/100.4°F), please use medical grade thermometer to verify.
- This product is not designed for use in medical evaluations. The forehead temperature is served as a reference only. It cannot be a judgment on fever.
- The use to this thermometer is not intended to be a substitution for consultation with your physician.
- To receive more accurate reading, the recommend measuring distance to forehead is 2~5 cm (0.8~2.0"). Please keep the forehead dry and clean to ensure the accuracy.
- For the best results, slowly scan from forehead towards the ear. Often the area just above the temple provides the best measurement because of it's proximity to the temporal artery.

- According to the clinical study, the reading of forehead measurement delivered by this product is equivalent to the oral temperature.
- Please be aware of the normal range of human body temperature varies due to an individual metabolism rate, the time of day, muscular activity, food intake or other factors.
- When serve the meter in ambient temperature 25±2°C (77±4°F), the meter will deliver the reading with the accuracy as stated. But when the ambient temperature changes drastically for example, (from outdoor to indoor), please leave it 30 minutes and start the service.
- The area temperature differs from the internal forehead temperature. To obtain the internal temperature always use the Forehead mode.

Please make sure to select the Forehead mode for an internal temperature reading and the Household mode for an external area reading(bottle, bath, room..).

2-2 Safety symbols

 Dangerous, refer to this manual before using the meter.

 CE Certification.

RoHS Restrict to use of six substances within electrical and electronic equipment (EEE), thereby contributing to the protection of human health and the environment.

 The device may not be disposed of with the trash. It promotes the re-use recycling and other forms of recovery of used materials and components, and to improve the environmental performance of all operators (manufacturers, traders,treatment facilities) involved in the life cycle of products. Dispose of the product appropriately in accordance with the regulations in force in your country.

3.Specifications


	Forehead	Household
Temperature Range	35.0~42.0°C (95.0~107.6°F)	0~100.0 °C (32.0~212.0°F)
Operation Temp.	25.0+2.0°C (77.0+4.0°F)	10.0~40.0°C (50.0~104.0°F)
Accuracy	+0.5°C (+1.0°F)	+1.0°C (+2.0 °F)
Distance/Spot Ratio	2~5cm (0.8~2.0")	8:1
Response Time	1 sec.	0.5 sec.
Thermopile	5~14 μm	
Repeatability	±1°C (±2°F)	
Resolution	0.1°C (0.1°F)	
Auto Power Off	Automatically after approx 20 sec.	
Storage Temp.	-10~60°C (14~140°F)	
°C / °F Switchable	YES	
Backlight	YES	
Battery Type	9V(006P, IEC6F22, NEDA1604)	
Battery Life	16hrs	
Dimension	156 X 100 X 32mm (6.1" x4.0"x1.3")	
Weight	135g Approx. (4.76 oz.)	
Accessory	9V Battery Instruction manual	

4.Operation of Instrument

4-1 Quick Start

To measure a temperature, point the unit at the target you want to measure, pull the trigger and hold. Be sure to consider the target area inside the angle of vision of this instrument.

The unit is powered by 9V battery and displays temperatures in either °C or °F. The user has to replace the battery when the battery voltage drops below the voltage for reliable operation and at the same time the low battery symbol  will appear.

5.Maintenance

Cleaning the lens: Blow off loose particles using clean compressed air. Gently brush remaining debris away with a camel's hair brush. Carefully wipe the surface with a moist cotton swab. The swab may be moistened with water.

NOTE:

DO NOT use solvents to clean the glass lens.

Cleaning the housing:

Use soap and water on a damp sponge or soft cloth.